

## ABSTRACT

A device for removing earwax is provided. With this device, it is not necessary to replace an earwax guard tip, and a user can remove earwax from the earwax guard tip easily on his/her own

5 The device is comprised of a case 1, a brush 2 accommodated in the case 1, a case 3, and a liquid retaining member 4 for absorbing disinfecting alcohol accommodated in the case 3. A user applies the earwax guard tip C to the liquid retaining member 4, lets the alcohol soak into the earwax, and thereafter rubs the earwax guard tip C against the brush 2 while moving the in-the-ear hearing aid H back and forth, left and right, for example.